

A HIGH-THROUGHPUT BACULOVIRUS PIPELINE ENABLING STRUCTURAL DETERMINATION OF INTEGRAL MEMBRANE PROTEINS

Nicola Burgess-Brown
PI Biotechnology

SGC, Nuffield Department of Clinical Medicine
University of Oxford, United Kingdom

Miptec, 24 – 25 September 2012




SGC Toronto





SGC Oxford


SGC PIPELINE



BACULOVIRUS EXPRESSION PROCESS




- 50-100 ml expression
- IMP pellets used for 12-detergent screening step
- Soluble protein pellets are discarded
- These could be purified for protein characterisation and/or crystallisation



METHOD DEVELOPMENT – PHYNEXUS PHYTIPS

PhyNexus Robot and Tip



Hopkins *et al* (2010) *J. Struct. Biol.*


- Received on loan for 8 months
- Chromatography using automated bench top platforms (for up to 12 samples in parallel)
- Selection of column bed sizes and chemistries (affinity, ion exchange, gel filtration, normal phase and reverse phase)
- Enhanced capture and elution of target proteins (in small volumes = higher conc.)

Our aims:


- To screen multiple constructs in crystallisation trials or biophysical characterisation in parallel from small volumes of culture
- To make use of the 50–100 ml pellets resulting from baculovirus amplifications
- To automate the purification process and free up time for other work
- To increase throughput in integral membrane protein (IMP) detergent screening

TESTING PHYTIPS ON EXISTING TARGETS


- Initial tests in *E. coli* on previously solved proteins
- Methods: Sonication (+ Benzonase), IMAC purification (eluted in 3x RBV) followed by desalting on spin columns
- 50 ml volume was the minimal quantity required for crystallisation experiments
- Low expressing proteins were not purified sufficiently
- Med/high expressing proteins gave concentrations up to 20 mg/ml




PHYTIPS FOR NOVEL TARGETS (*E. COLI*)



- 8 constructs of a novel myosin protein tested in parallel
- In less than a week, 8 constructs purified from 50 ml cultures and **2** of these resulted in crystals



- 8 constructs of a novel demethylase tested in parallel
- 3 successfully purified, but only **1 construct** resulted in crystal hits (conditions similar to those used for other demethylases) – diffracted to 2.6 Å at Diamond




- A novel 3 protein kinase complex purified to 8.3 mg/ml gave crystal hits

PHYTIPS FOR SCREENING IMPs – PRELIMINARY DATA

- **Aim:** To test the use of PhyTips for automated purification of IMPs for monodispersity studies by gel filtration
- **Method:** 100–250 ml insect cell culture tested, membrane prep, 20 μ l bed volume of Ni-NTA resin or Talon resin

ZMPSTE24 Purified Manually and using PhyTips



Conclusion:

- The purity resulting from PhyTips was comparable to the manual batch method
- Total yield recovered from the PhyTips was 2-fold higher for an equivalent resin bed volume (RBV)

Conclusions so far:

- Buffer screen for 3 ABC transporters gave useful conditions and freed up researcher time
- Also suitable for larger lysate volumes (e.g. 10 ml purified in 3 h)




Manual drip columns

PHYNEXUS TIPS – CONCLUSIONS SO FAR

- Potential for higher throughput expression screening to crystal trials (using robotics)
- Allows us to obtain small amounts of protein **at high concentration**, because of the small elution volume
- Suitable for med/high expressers from 50 ml bacterial expression
- Useful for screening multiple constructs of 1 target in parallel
- Potential for increasing throughput of IMP screening (buffers, detergents, ligands, etc)

Still to test/optimize:

1. Baculovirus pellets
2. On column tag cleavage
3. HTP desalting
4. Other resins/chemistries
5. Using our Hamilton robot with PhyTips



ACKNOWLEDGEMENTS



Biotechnology

Nicola Burgess-Brown

Georgina Berridge

Rod Chalk

Leela Shrestha

Solenne Goubin

Claire Strain-Damerell

Grazyna Kochan

Pravin Mahajan

Shubhashish Mukhopadhyay

Emilia Hapka

Aleksandra Szykowska

Integral Membrane Protein Group

Liz Carpenter

Ashley Pike

Andrew Quigley

Liang Dong

Yin Dong

Xiaobing Xia

Chitra Shintre

Mariana Grieben

Everyone else at Oxford SGC

Ex-SGC members

Juin-In Kim

Sanjan Das

Beatrice Li

Jing Yang

Ivan Alfano

James Bray

Binesh Shrestha

Alastair Barr

Chas Bountra (CSO Oxford)

Aled Edwards (CEO)

Cheryl Arrowsmith (CSO Toronto)

Opher Gileadi (Genome Integrity & previous Biotech PI)

Stefan Knapp (Kinases, Epigenetics)

Alex Bullock (Growth Factor Signalling)

Paul Brennan (Medicinal Chemistry)

Frank von Delft (Crystallography)

Brian Marsden (Research Informatics)

Wen Hwa Lee (Scientific Coordinator)

Susanne Müller-Knapp (Epigenetics)

Wyatt Yue (Organelle Biogenesis)

Scott Thomson (Administrator)

Udo Oppermann (Metabolic Enzymes, Epigenetics)

Johan Weigelt (former CSO Stockholm)

Diamond Light Source Ltd

FUNDING PARTNERS

The Canadian Institutes for Health Research, the Canada Foundation for Innovation, Genome Canada, GlaxoSmithKline, Lilly Canada, the Novartis Research Foundation, Pfizer, Takeda, the Ontario Ministry of Economic Development and Innovation, and the Wellcome Trust.

CONTACT: nicola.burgess-brown@sgc.ox.ac.uk

ABOUT THE SGC

The SGC is an international public-private partnership (UK charity number 1097737) that aims to carry out basic science of relevance to drug discovery, placing all information, reagents and know-how into the public domain without restriction.