

PhyNexus AutoPlasmid MMG Kits

To reorder, call (408) 267-7214 or email: sales@phynexus.com

Maxiprep Kits (16 samples)		Part Number: DPM 16-03-70-01
Description		QTY
Tips	PhyTip Columns (3 mL resin)	16 x
	Pipette Transfer Tips (20 mL)	48 x
Buffers	Resuspension Buffer (48 mL in 60 mL bottles)	16 x
	Lysis Buffer (57 mL in 90 mL bottles)	16 x
	Precipitation Buffer (68 mL in 90 mL bottles)	16 x
	Wash Buffer B (130mL in 500 mL bottle; user adds 300 mL EtOH for a total of 430 mL)	1 x
	Wash Buffer (1184 mL in 1 gallon bottle; user adds 2200mL EtOH for a total of 3384 mL)	1 x
Parts	Elution Buffer (20 mL)	16 x
	Elution tubes (5 mL snap cap)	16 x
	Sample plates	16 x
	RNase A tubes (20 mg)	16 x

Megaprep Kits (16 samples)		Part Number: DPM 16-06-70-01
Description		QTY
Tips	PhyTip Columns (6 mL resin)	16 x
	Pipette Transfer Tips (20 mL)	48 x
Buffers	Resuspension Buffer (144 mL in 250 mL bottles)	16 x
	Lysis Buffer (171 mL in 250 mL bottles)	16 x
	Precipitation Buffer (204 mL in 250 mL bottles)	16 x
	Wash Buffer B (130mL in 500 mL bottle; user adds 300 mL EtOH for a total of 430 mL)	2 x
	Wash Buffer (1184 mL in 1 gallon bottle; user adds 2200mL EtOH for a total of 3384 mL)	2 x
Parts	Elution Buffer (20 mL)	16 x
	Elution tubes (5 mL snap cap)	16 x
	Sample plates	16 x
	RNase A tubes (60 mg)	16 x

Gigaprep Kits (8 samples)		Part Number: DPM 08-06-70-01
Description		QTY
Tips	PhyTip Columns (6 mL resin)	16 x
	Pipette Transfer Tips (20 mL)	48 x
Buffers	Resuspension Buffer (144 mL in 250 mL bottles)	16 x
	Lysis Buffer (171 mL in 250 mL bottles)	16 x
	Precipitation Buffer (204 mL in 250 mL bottles)	16 x
	Wash Buffer B (130mL in 500 mL bottle; user adds 300 mL EtOH for a total of 430 mL)	2 x
	Wash Buffer (1184 mL in 1 gallon bottle; user adds 2200mL EtOH for a total of 3384 mL)	2 x
Parts	Elution Buffer (20 mL)	16 x
	Elution tubes (5 mL snap cap)	16 x
	Sample plates	16 x
	RNase A tubes (60 mg)	16 x