



PhyNexus, Inc.

IMAC PHYTIP COLUMNS MSDS

MATERIAL SAFETY DATA SHEET

PhyNexus, Inc.
3670 Charter Park Drive, Suite A
San Jose, CA 95136
United States

Preparation Date: Jun 21, 2016
Printed Date: Jun 21, 2016
Technical Phone: 408-267-7214

SECTION 1. CHEMICAL IDENTIFICATION

NAME: IMAC PHYTIP COLUMNS
SUBSTANCE: POLYMERS WITH PORTION COATED WITH GLYCEROL
CAS#: None

SECTION 2. COMPOSITION/INFORMATION ON INGREDIENTS

CAS#: 9003-07-0 POLYPROPYLENE
MF: N.A.
EC Number: Not Assigned
SYNONYMS: N.A.
Polypropylene is considered non hazardous.

CAS#: 9012-36-6 POLYSACCHARIDE AGAROSE
MF: N.A.
EC Number: Not Assigned
SYNONYMS: Agarose
Polysaccharide agarose is considered non hazardous

CAS#: 56-81-5 GLYCEROL 10%
Chemical Formula: $C_3H_5(OH)_3$
SYNONYMS: 1,2,3-PROPANETRIOL; GLYCERIN; GLYCOL ALCOHOL; Glycerol is generally regarded as a safe or slightly hazardous material. It is slightly irritating and slightly toxic. It is flammable.

CAS#: 7732-18-5 WATER 10%

SECTION 3. HAZARDS IDENTIFICATION

NFPA Ratings (SCALE 0-4): HEALTH=0-1 FIRE=0-1 REACTIVITY=0
NON HAZARDOUS OR SLIGHTLY HAZARDOUS
CARCINOGEN STATUS:
-OSHA: No.
-NTP: No.
-IARC: No.



IMAC PHYTIP COLUMNS MSDS

SECTION 4. FIRST-AID MEASURES

SEEK MEDICAL ASSISTANCE FOR ANY OVEREXPOSURE.

SKIN: Wash thoroughly with soap and water.

EYES: Immediately flush thoroughly with water for at least 15 minutes.

INGESTION: If conscious drink copious amounts of water. Never give anything by mouth to an unconscious person. Obtain medical attention.

SECTION 5. FIRE FIGHTING MEASURES

EXTINGUISHING MEDIA: Use extinguishing agents appropriate for surrounding fire.

SPECIAL FIREFIGHTING PROCEDURES: Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

UNUSUAL FIRE AND EXPLOSION HAZARDS: None, not a fire hazard.

SECTION 6. ACCIDENTAL RELEASE MEASURES

Absorb material and collect in an appropriate container for disposal. Wash spill site with soapy water.

SECTION 7. HANDLING AND STORAGE

Store and handle in accordance with all current regulations and standards. Keep separated from incompatible substances. Store between 4-20°C, do not freeze.

SECTION 8. EXPOSURE CONTROLS AND PERSONAL PROTECTION

EXPOSURE LIMITS: No occupational exposure limits established.

CHEMICAL SAFETY GOGGLES.

SAFETY SHOWER AND EYE BATH.

COMPATIBLE CHEMICAL-RESISTANT GLOVES (PVC, RUBBER OR NITRILE).

AVOID CONTACT WITH EYES, SKIN AND CLOTHING.

KEEP TIGHTLY CLOSED.

STORE BETWEEN 4-20°C.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

SPECIFIC GRAVITY: N.A.

APPEARANCE: Polymer and Colorless liquid.

ODOR: Odorless.

BOILING POINT: N.A.

WATER SOLUBILITY: Portions are Water soluble. Polymers are insoluble.



IMAC PHYTIP COLUMNS MSDS

SECTION 10. STABILITY AND REACTIVITY

REACTIVITY: Stable at normal temperatures and pressure.
INCOMPATIBILITIES: Strong oxidizing agents.
HAZARDOUS POLYMERIZATION: Will not occur.
DECOMPOSITION PRODUCTS: Carbon Dioxide, Carbon Monoxide.

SECTION 11. TOXICOLOGICAL INFORMATION

INHALATION: No Data Available.
SKIN CONTACT: No Data Available.
EYE CONTACT: No Data Available.
INGESTION: No Data Available.

SECTION 12. ECOLOGICAL INFORMATION

Data Not Yet Available.

SECTION 13. DISPOSAL CONSIDERATIONS

Dispose in accordance with all applicable regulations.

SECTION 14. TRANSPORT INFORMATION

No classification assigned (DOT, Canadian, Land, Air & Maritime).

SECTION 15. REGULATORY INFORMATION

TSCA INVENTORY STATUS: No.
TSCA 12 (b) EXPORT NOTIFICATION: Not listed.
CERCLA SECTION 103 (40CFR302.4): No.
SARA SECTION 302 (40CFR355.30): No.
SARA SECTION 304 (40CFR355.40): No.
SARA SECTION 313 (40CFR372.65): No.
SARA HAZARD CATEGORIES, SARA SECTIONS 311/312 (40CFR370.21):
 ACUTE: No.
 CHRONIC: No.
 FIRE: No.
 REACTIVE: No.
 SUDDEN RELEASE: No.
OSHA PROCESS SAFETY (29CFR1910.119): No.
STATE REGULATIONS:
 CALIFORNIA PROPOSITION 65: No.
CANADIAN REGULATIONS: Not assigned.
EUROPEAN REGULATIONS: Not assigned.



PhyNexus, Inc.

IMAC PHYTIP COLUMNS MSDS

SECTION 16. ADDITIONAL INFORMATION

The information provided herein is believed to be accurate at the time of preparation or prepared from sources believed to be reliable. Users have the responsibility to comply with all laws and procedures applicable to the safe handling and use of the product and to determine the suitability of the product for its intended use. PHYNEXUS MAKES NO WARRANTY, EXPRESSED OR IMPLIED, CONCERNING THE PRODUCT OR THE MERCHANTABILITY OR FITNESS THEREOF FOR ANY PURPOSE OR CONCERNING THE ACCURACY OF ANY INFORMATION PROVIDED, AND IN NO WAY SHALL PHYNEXUS BE LIABLE FOR ANY CLAIMS, LOSSES, OR DAMAGES OF ANY THIRD PARTY OR FOR LOST PROFITS OR ANY SPECIAL, INDIRECT, INCIDENTAL, CONSEQUENTIAL OR EXEMPLARY DAMAGES HOWSOEVER ARISING.

PhyNexus, Inc.

3670 Charter Park Drive, Suite A, San Jose, CA 95136

Tel: 408.267.7214 Fax: 408.267.7346

www.phynexus.com